**TeDan Surgical
Innovations**

Phantom UL Anterior to the Psoas (ATP) Surgical Access System

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TABLE OF CONTENTS

1. INTRODUCTION	3
1.1 PRODUCT DESCRIPTION	3
1.2 SURGICAL SETTING	3
1.3 INDICATIONS FOR USE.....	3
1.4 INTENDED PATIENT GROUPS.....	3
1.5 INTENDED USER GROUPS.....	3
1.6 CONTRAINDICATION	3
1.7 INTENDED PURPOSE	3
1.8 INTENDED CLINICAL BENEFITS.....	3
1.9 SYMBOLS GLOSSARY	3
1.10 PRIOR TO USE	3
2. SET UP.....	3
2.1 PATIENT POSITIONING	3
2.2 FLUOROSCOPY	3
2.3 INCISION PLANNING	4
2.4 ARM ASSEMBLY	4
2.5 RETRACTOR RING ASSEMBLY - CHOOSE 3/4 OR HALF RING	4
2.5.1 3/4 RING	4
2.5.2 HALF RING	5
2.6 RETRACTOR RING TO ARTICULATING ARMS ASSEMBLY	5
2.7 BLADE SELECTION	5
2.8 BLADE TO BLADE HOLDER ASSEMBLY	6
2.9 RING CLAMP TO RETRACTOR RING ASSEMBLY AND POSITIONING	6
2.10 ATP SCREWS DEPLOYMENT	6
3. LIGHTING SYSTEM ASSEMBLY - CHOOSE SHIM TIP OR PRONGED BUSHING	7
3.1 SHIM TIP (ML-0068)	7
3.2 PRONGED BUSHING (ML-0058)	7
4. OPTIONAL: L5-S1 ASSEMBLY	7
4.1 OPTIONAL: FLEX ARM ASSEMBLY	7
4.2 OPTIONAL: FLEX ARM TO ARTICULATING ARM ASSEMBLY	8
4.3 OPTIONAL: RETRACTOR BLADE TO FLEX ARM ASSEMBLY	8
5. DISASSEMBLY	8
6. CLEANING AND STERILIZATION WARNINGS	8
6.1 REPROCESSING INSTRUCTIONS	8
7. CLEANING INSTRUCTIONS	8
7.1 POINT OF USE	8
7.2 SUPPLIED AND EQUIPMENT NEEDED FOR CLEANING	9
8. CLEANING	9
8.1 MANUAL CLEANING INSTRUCTIONS	9
8.2 AUTOMATED CLEANING INSTRUCTIONS	10
9. POST CLEANING INSPECTION	10
9.1 SURGICAL ARM INSPECTION	10
10. STERILIZATION INSTRUCTIONS	11
10.1 LUBRICATE	11
10.2 STERILIZE	11
11. STORE	11
12. SERIOUS INCIDENTS REPORTING	11
13. PRODUCT LIFETIME	11
14. WARNINGS/PRECAUTIONS	12

1. INTRODUCTION:

1.1 PRODUCT DESCRIPTION:

The TeDan Surgical Innovations, Inc. (TSI) Phantom UL™ Anterior to the Psoas (ATP) Surgical Access System enables the surgeon to access and work on the lumbar disc space between L1 and L5 as well as L5 and S1; to remove the disc, prepare end plates, and insert and fixate an implant.

1.2 SURGICAL SETTING:

Always ensure use of TSI instruments in a strictly sterile clinical environment.

1.3 INDICATIONS FOR USE:

Anterior to the Psoas surgery for lumbar levels L1 to S1.

1.4 INTENDED PATIENT GROUPS:

Patients for which lumbar spine surgery from an anterior to the psoas approach is required.

1.5 INTENDED USER GROUPS:

Licensed practicing physicians who have been informed on the use of the ATP Surgical Access System. It is also expected that clinical assistants are knowledgeable in setting up the ATP Surgical Access System.

1.6 CONTRAINDICATION:

No known contraindications exist.

1.7 INTENDED PURPOSE:

The ATP Surgical Access System is intended to facilitate minimally invasive access to the lumbar spine from an oblique angle providing retraction of relevant anatomy out of the surgical area.

1.8 INTENDED CLINICAL BENEFITS:

The ATP Surgical Access System enables the surgeon to access and work on the disc space; to remove the disc, prepare end plates, and insert and fixate an implant.

1.9 SYMBOLS GLOSSARY:

See www.tedansurgical.com/symbolsglossary for the Symbols & Device Markings Glossary.

1.10 PRIOR TO USE:

- Clean and sterilize all ATP components according to the Cleaning and Sterilization Instructions provided below.
- For instruments with moving parts, lubricate joints with a steam permeable, water soluble instrument lubricant prior to sterilization.
- Articulating arm inspection before use: Inspect entire assembly for damage. Hold arm assembly and turn central tightening knob clockwise. Check to make sure that arm is rigid at all three joints.

2. SET UP:

2.1 PATIENT POSITIONING:

- Place the patient in the lateral decubitus (left side up) position (Figure 1).

2.2 FLUOROSCOPY:

- Perform fluoroscopy to identify affected disc.


Figure 1

2.3 INCISION PLANNING:

For L1-L5	For L5 - S1
- Based on fluoroscopy images taken, mark the end plate of affected discs to identify where the incision will be made. - Use markings to determine ring placement in steps below. The ring should be centered around the incision.	- Localize the affected disc space in the L5 - S1 region and mark accordingly. - Use markings to determine ring placement in steps below. Ring should be centered or slightly cephalad in as close to the level of the umbilicus.

2.4 ARM ASSEMBLY:

- Hand thread the quick connection of the **Articulating Arm** (UL-0110) onto the Articulating Arm column (This is more easily done by holding the post vertical, thread down, with the articulating portion flat on the table). The post will guide itself into place.
- Use the **Arm Wrench** (HS-0110) to tighten the quick connection so that it is secure. (Figure 2) Repeat the above steps for the second arm.
- Hold assembled Articulating Arms at desired place on the OR table to establish placement for the **Table Clamps** (ML-0021).
- Suggested positions for setup include a clamp positioned cephalic, left posterior and one clamp positioned caudally, right anterior. (Figure 3)
- Loosen knob of **Table Clamp** (ML-0021). Attach **Table Clamp** to the surgical rail over the sterile drape in the established position based on the previous step.
- Tighten knob to secure **Table Clamp** at the desired location. Using the same method, attach the second **Table Clamp** to the opposite side of the table.
Note: Table clamps should be deployed so that one table clamp is positioned as far caudally with the other as far cephalic as possible while still allowing the arm clamp to attach to the retractor ring.
- Insert an **Articulating Arm** (UL-0110) into each of the **Table Clamps**. Position the **Articulating Arms** and **Table Clamps** so that they are angled as low as possible and will not impede the surgeon's movement. (Figure 4)

Note: When loosening, do not force the knob of the articulating arm past stop. Doing so could damage the ball joint and affect the rigidity of the articulating arm.


Figure 2

Figure 3


Figure 4

2.5 RETRACTOR RING ASSEMBLY- CHOOSE 3/4 OR HALF RING:

2.5.1 ¾ RING :

- Clear enough space on a flat surface for assembly of the **ATP Hinged ¾ Ring** (UL-0100). Each of the following **Retractor Ring** steps must be conducted on a flat surface to ensure that the screws and forks are properly aligned.
- ATP Hinged ¾ Ring** consists of two ring segments. To loosen the screw, turn the screw counterclockwise no more than two times on the ¾ ring segment using the **Hex Tool** (ML-0505).
- Using the **Hex Wrench** (AL-0106) or by hand, loosen the hex screw of 1/4 ring segment (Turn the hinge hex counterclockwise no more than two times on each ring segment). Assemble the **ATP Hinged ¾ Ring** by pushing the two segments together so that the ¾ ring's screw side interlocks with the fork side of the half ringed portion.

Note: Connection point of the ATP Hinged 3/4 Ring should be brought to the edge of flat surface to provide clearance for Hex Tool during assembly.

- With the **Hex Tool** (ML-0505), tighten the screw onto the half ringed portion to secure the assembly. (Figure 5 and 6)

Note: To prevent damage to the screws and forks on the connection points of the ATP Hinged 3/4 Ring, pull segments directly apart. Do not cant or twist the rings. Do not force screws past stop.


Figure 5


Figure 6

2.5.2 HALF RING:

- Clear enough space on a flat surface for assembly of the **ATP Half Ring**. Each of the following **Retractor Ring** steps must be conducted on a flat surface to ensure that the screws and forks are properly aligned.
- **ATP Hinged ¾ Ring** (UL-0100) consists of two ring segments. For the **ATP Half Ring** set up, use only the half ring portion of the **ATP Hinged ¾ Ring**.
- One end of the half ring portion has a fork side. Assemble the **ATP Half Ring** by pushing the fork side together with the **Ring End Cover, Fork Side** (UL-0102) until it interlocks with the fork side of the half ringed portion. To loosen the screw, turn the screw counterclockwise no more than two times on the **Ring End Cover, Fork Side** using the **Hex Tool** (ML-0505).

Note: To prevent damage to the screws and forks on the connection points of the ATP Half Ring, pull segments directly apart. Do not cant or twist the rings. Do not force screws past stop.

- With the **Hex Tool** (ML-0505), tighten the screw onto the half ringed portion to secure the assembly.
Note: Connection point of the ATP Half Ring should be brought to the edge of flat surface to provide clearance for Hex Tool during assembly.

2.6 RETRACTOR RING TO ARTICULATING ARMS ASSEMBLY:

- Using the **Hex Wrench** (AL-0106), unscrew the clamps of both **Articulating Arms**. Insert the **Retractor Ring** into the clamps and tighten the clamps around the ring where marked "Attach Rigid Arm Here". (Figure 7)

Note: Choose the appropriate side of the Retractor Ring, per the Level Markings on either side of the Retractor Ring (L1-L5 or L5/S1)

Note: Ensure the hex wrench and screw are aligned with the threaded hole.

Attempting to insert the screw at an angle may cause cross-threading. To avoid damage to mechanism, do not force screw mechanism past stop.

- Adjust the **Articulating Arms** to position the **Retractor Ring** in accordance with the markings that indicate the Caudal and Cephalic direction. To prevent movement of the **Retractor Ring**, ensure that the knobs of the **Articulating Arms** are completely tightened. Using the Hex Wrench (AL-0106) or by hand, tighten the hex screw on the 1/4 portion of the hinged ring.

Note: The 1/4 portion of the hinged ring allows the retractor ring to bend and be adjusted to conform to the patient's anatomy.

Figure 7

2.7 BLADE SELECTION

For L1-L5	For L5-S1
• Cannulated ATP Blade: Use two of these blades to retract the psoas posteriorly. (Figure 8)  Figure 8 • Nose Tip ATP Blade: Use one blade to retract/position the vessels on the contralateral side of the disc space. Position the blade along the disc space, and then slide between the disc and vessels to retract contralaterally from the vertebral body. (Figure 9)  Figure 9	• Cannulated ATP Blade: Use one of these blades under the iliac vein. This step protects and lateralizes the left common iliac vein and artery. A second cannulated blade can be used at the bifurcation at L5, this blade may be pinned to the L5 vertebral body. (Figure 8) • Nose Tip ATP Blade: Use one blade to retract/position the vessels on the contralateral side of the disc space. Position the blade along the disc space, and then slide between the disc and vessels to retract contralaterally from the vertebral body. (Figure 9) • ALIF Blade: Use one blade to retract the renal structures in the sacral direction. (Figure 10)  Figure 10

2.8 BLADE TO BLADE HOLDER ASSEMBLY:

- **Retractor Blades** have a “dual fixation” connection which allows the blade to remain fixed or to rotate. To attach the **Retractor Blade** to the **Blade Holder** (such as AL-0105), insert the **Retractor Blade** through the **Blade Holder** connection so that the first pin in the connector can pass through the groove in the **Blade Holder**. (Figure 11)
- When the first pin is engaged with the **Blade Holder**, the **Retractor Blade** is attached and able to rotate.
- If a fixed connection is preferred, the second pin should be aligned with a groove, and the **Retractor Blade** pushed upward. When the second pin is engaged in the **Blade Holder**, the **Retractor Blade** will remain fixed.


Figure 11

2.9 RING CLAMP TO RETRACTOR RING ASSEMBLY AND POSITIONING:

For L1-L5	For L5-S1
• Suggested placement of blades in the incision: Cannulated ATP Blades (A), Nose Tip ATP Blade (B) (Figure 12) Figure 12	• Suggested placement of blades in the incision: Nose Tip ATP Blade (C), Cannulated ATP Blades (D), ALIF Blade (E) (optional, only needed for four-way retraction in the sacral direction) (Figure 13) Figure 13

- **Ring Clamps** (AL-0103) attach to any point on the ring. To attach, loosen knob of the **Ring Clamp** and press the larger of the two openings onto the **Retractor Ring**; the **Ring Clamps** will snap onto the **Retractor Ring**. Press the shaft of the **Blade Holder** into the smaller opening of the **Ring Clamp** to attach the two components. The **Ring Clamp/Blade Holder** Assembly can be freely adjusted until the knob of the clamp is tightened. Additional **Ring Clamps** and **Blade Holders** may be added to the **Retractor Ring** as desired. (Figure 14)
- To angulate the **Retractor Blades**, use the **Hex Wrench** (AL-0106) to loosen the swivel mechanism on the **Blade Holder** and manipulate the **Retractor Blade** by hand to the desired position. To secure the **Retractor Blade**, tighten the swivel mechanism with the **Hex Wrench**. (Figure 15)

Note: It is only necessary to rotate the mechanism counterclockwise until the swivel mechanism moves freely. To avoid damage to mechanism, do not force screw mechanism past stop. Do NOT use Hex Wrench to pivot the blade.


Figure 14


Figure 15

2.10 ATP SCREWS DEPLOYMENT:

- Choose **ATP Screws** (UL-0360 to UL-0378) based on blade length, use **Lateral Screw Driver** (such as ML-0055) to screw into vertebral body. Screw should go into a cannulation on the **Cannulated ATP blade**. (Figure 16)
- Only one **ATP Screw** should be used per **Cannulated ATP blade**.


Figure 16

For L1-L5	For L5-S1
After each Cannulated ATP Blade is placed in the desired position, ATP Screws should be added into the desired cannulations for the Cannulated ATP Blades.	Each ATP Screw should be placed after each Cannulated ATP Blade is placed in the desired position. One ATP Screw should be inserted into the S1 vertebral body. The second ATP Screw should be inserted into the L5 vertebral body.

3. LIGHTING SYSTEM ASSEMBLY- CHOOSE SHIM TIP OR PRONGED BUSHING:

3.1 SHIM TIP (ML-0068):

- Plug the end of the **Light Cable, Shim Tip, 10', Non-sterile (ML-0068)** into the **LED Light Source (ML-0051)**
- Once the **Light Cable** is inserted into the **LED Light Source**, open the light source shutter to confirm that the **Light Cable** is transmitting light. Once confirmed, close the light source shutter. The light source shutter must remain closed until the **Light Cable** is inserted into the **Retractor Blades**.
- Slide the shim tips of the **Light Cable, Shim Tip, 10', Non-sterile (ML-0068)** into the shim channels of the **Retractor Blades**. (Figure 17)

Note: Light cables have a malleable section that allows the user to fix each tip in place.

- Open the **LED Light Source** shutter and confirm that the **Light Cable** is transmitting light.
- The shutter on the **LED Light Source** may be adjusted to provide the amount of light preferred.
- Close the **LED Light Source** shutter prior to removing the **Light Cable** from the **Retractor Blades**.


Figure 17

3.2 PRONGED BUSHING (ML-0058):

- Plug the end of the **Secure Light Cable, Pronged Bushing, 10', Non-sterile (ML-0058)** into the **LED Light Source (ML-0051)**
- Once the **Light Cable** is inserted into the **LED Light Source**, open the light source shutter to confirm that the **Light Cable** is transmitting light. Once confirmed, close the light source shutter. The light source shutter must remain closed until the **Light Cable** is inserted into the **Retractor Blades**.
- Gently insert the bifurcated ends into the holes located on top of the **Retractor Blades**. (Figure 18)

Note: Light cables have a malleable section that allows the user to fix each tip in place.

- Open the **LED Light Source** shutter and confirm that the **Light Cable** is transmitting light.
- The shutter on the **LED Light Source** may be adjusted to provide the amount of light preferred.
- Close the **LED Light Source** shutter prior to removing the **Light Cable** from the **Retractor Blades**.


Figure 18

4. OPTIONAL L5-S1 ASSEMBLY:

- Turn the flex arm tightening knob clockwise and ensure the arm is sufficiently rigid for intended use.
- When loosened, check cable between links for fraying. Normal use will eventually cause wear to the steel tensioning cable. If cable shows any frayed or broken wires, the flex arm needs to be replaced immediately

4.1 OPTIONAL: FLEX ARM ASSEMBLY:

- To attach the **Quick Connect Blade Holder** (such as UL-0107) to **Flex Arm** (such as UL-0106), ensure the **Flex Arm** collar is in its unlocked position. (Figure 19)
- Fully insert the **Quick Connect Blade Holder** shaft into the collar of the **Flex Arm**.
- Twist the collar clockwise to lock the **Quick Connect Blade Holder** into place.
- Check to see that the **Quick Connect Blade Holder** is secure and the collar is in the locked position.


Figure 19

4.2 OPTIONAL: FLEX ARM TO ARTICULATING ARM ASSEMBLY:

- Loosen the clamp using the small purple knob of the **Flex Arm** (such as UL-0106).
 - Position the jaws around the segments of either **Articulating Arm** (UL-0110).
 - Tighten small purple knob to affix the **Flex Arm** to the **Articulating Arm**. (Figure 20)
- Note: Do not force past stop.*


Figure 20

4.3 OPTIONAL: RETRACTOR BLADE TO FLEX ARM ASSEMBLY:

- ALIF Blades** have a “dual fixation” connection which allows the **ALIF Blade** to remain fixed or to rotate. To attach the **ALIF Blade** to the **Quick Connect Blade Holder** (such as UL-0107), insert the **ALIF Blade** through the **Quick Connect Blade Holder** connection so that the first pin in the connector can pass through the groove in the **Quick Connect Blade Holder**.
- When the first pin is engaged with the **Quick Connect Blade Holder**, the **ALIF Blade** is attached and able to rotate.
- If a fixed connection is preferred, the second pin should be aligned with one of the grooves, and the **ALIF Blade** pushed upward. When the second pin is engaged in the **Quick Connect Blade Holder**, the **ALIF Blade** will remain fixed. (Figure 21)
- Tighten large purple knob to fixate positioning of the **Flex Arm**.


Figure 21

5. DISASSEMBLY:

- Completely loosen the **Articulating Arm** clamp screws, do not over-loosen or force the screw past stop.
 - Disassemble the **Retractor Ring** by turning the screws counterclockwise twice with the **Hex Tool**.
 - It is important that the **Retractor Ring** is disassembled by pulling the segments directly apart on a flat surface.
- Note: To prevent damage to the screws and forks on the connection points of the Retractor Ring, pull rings directly apart. Do not cant or twist the rings. Do not over-loosen or force the screws past stop.*
- Follow the above assembly steps (2 - 4.3) in reverse order to properly disassemble the system.
 - ATP screws are single use only products. After the procedure, dispose of screws in an approved sharps container.

6. CLEANING AND STERILIZATION WARNINGS:

Please note:

- The color of TSI's Titanium and Aluminum instruments may vary due to the anodizing process or alloy used. Shading or loss of color may also occur after sterilization. This is not a defect in the instrument or material and will not affect the performance of your high quality TSI instrument.
- When loosening, do not force the knob of the articulating arm past the stop. Doing so could damage the ball joint and affect the rigidity of the articulating arm (Figure 22).
- Forcing the pivoting mechanism past stop may cause damage to the device.

Figure 22

6.1 REPROCESSING INSTRUCTIONS:

The Cleaning and Sterilization Steps defined in this Instruction for Use, should be used for reprocessing of TSI products. Other methods used for reprocessing of TSI devices shall be validated by the user prior to implementation.

7. CLEANING INSTRUCTIONS:

7.1 POINT OF USE:

- Directly after use, remove coarse contamination from the instrument and keep the instrument moist for transit to the processing site. Prior to cleaning and sterilization do not use any fixing agents or hot water >104°F (>40°C) since this may lead to the fixation of residue and can interfere with the cleaning process.
- Where applicable, multi-component instruments should be disassembled for appropriate cleaning. Care should be exercised to avoid losing small components.
- Perform cleaning as soon as is reasonably practical following use. Allowing instrumentation to dry with fluids and debris may result in staining, corrosion, and increased difficulty in the removal of contaminants.
- During the transport of the instruments to the processing site, store contaminated instrument(s) securely in a closed container to avoid damage to the instrument and/or contamination of the environment.

7.2 SUPPLIES AND EQUIPMENT NEEDED FOR CLEANING:

- Ultrasonic Cleaner
- Pre-Cleaning and Manual Cleaning enzymatic, pH Neutral Cleaner: such as Prolystica 2x Concentrate Enzymatic Pre-Soak and Cleaner
 - Soft Nylon Bristle Brush 3-1000 Integra Miltec Premium Grade Nylon Bristle Brush should be used for exterior surface cleaning. For internal cleaning, a brush size appropriate to the lumen inner diameter (ID) should be used.
- Automated Cleaning enzymatic, pH Neutral Cleaner: such as Prolystica 2x Concentrate Enzymatic Pre-Soak and Cleaner
- 60 ml Syringe with Needle
- Soft Bristled Brush (ex. M16)
- Tap, Deionized (DI), Reverse Osmosis (RO), or filtered water for processing

8. CLEANING:

***Automated cleaning is not suitable for instruments with long lumens, ball joints, or braided cables (e.g. articulating surgical arms, and flexible surgical arms). Such instruments should only undergo manual cleaning.**

- For both Manual and Automated Cleaning, all instruments, except the Articulating Arm, should be cleaned in the open or unlocked position.
 - **Flexible Arms:** Turn tightening knob counterclockwise to loosen the internal cable to enable movement of the surgical arm beads to allow water or detergent to flow between each link prior to placing arm in sonicator.
 - **Articulating Arms:** Turn tightening knob clockwise to tighten ball joint prior to placing arm in sonicator. The knob should be tightened enough to prevent manipulation of the arm.
- An enzymatic, pH neutral cleaner is recommended.
 - Black coated and color anodized components may be negatively affected if aggressive cleaning mediums or appliances (e.g. extreme acidic/alkaline, abrasives) are used.
 - Exposure to chlorides or hydrogen peroxide may negatively affect the coating or colorization of the components.

For the purposes of this IFU, the table below defines cold, lukewarm, and hot water temperatures (per AAMI TIR12:2010).

Temperature Description	°Celsius	°Fahrenheit
Cold	<22°C	<72°F
Lukewarm	22°- 43°C	72°-110°F
Hot	>43°C	>110°F

8.1 MANUAL CLEANING INSTRUCTIONS:

1. Rinse each instrument individually with a steady stream of lukewarm tap water for a minimum of 2 minutes or until gross contaminants are removed.
 - For Flex Arms, slide the surgical arm beads to allow water to flow between each link.
2. Place each instrument into a sonicator containing a solution of enzymatic, pH neutral detergent prepared according to the manufacturer's instructions, using lukewarm tap water. Sonicate for 10 minutes.
 - Ensure lumens are flushed to remove air bubbles prior to sonication.
3. In a manual wash container, prepare an enzymatic, pH neutral detergent wash solution, per detergent manufacturer's instructions, using lukewarm tap water.
4. Transfer each instrument from the sonicator to the manual wash container, and fully immerse in the cleaning solution prepared in Step 3.
5. While still submerged, using a soft nylon bristled brush, such as M16, brush each instrument for a minimum of 1 minute, to remove any visible contamination and debris, paying particular attention to hard-to-clean areas such as crevices and joints. Repeat the process until all visible contamination is removed.
 - For blades with lumens, use appropriately sized lumen brush to scrub the interior of lumen and soft bristled brush, such as M16, to brush the exterior of the lumen. Using a 60mL syringe, flush lumens with 25mL of cleaning solution, a minimum of 3 times, or until no soil is visible.
 - For Flex Arms, turn the tightening knob counterclockwise to loosen the internal cable to enable movement of the surgical arm beads to allow cleaning solution fluid to flow between each link.

6. Rinse each instrument with lukewarm tap water for a minimum of 2 minutes, or until no visible soil remains.
 - For items with lumens, using a 60mL syringe with needle, flush the lumen with 25mL of lukewarm water 3 times.
 - For Flex Arms, slide the surgical arm beads to allow water to flow between each link.
7. Dry each instrument using clean, absorbent, lint-free wipes, or pressurized air, to remove excess rinse water.

8.2 AUTOMATED CLEANING INSTRUCTIONS:

Use only washer/disinfector machines that have been validated in accordance with ISO 15883.

Articulating Surgical Arms and Flexible Surgical Arms should not be cleaned using the Automated Method. These devices should be cleaned using the Manual Cleaning Process (detailed above).

1. Perform pre-cleaning to remove gross contaminants as follows:
 - a. Rinse with running, lukewarm, DI, RO, or filtered water for a minimum of 1 minute for each instrument to remove gross contaminants.
 - b. In a manual wash container, prepare an enzymatic, pH neutral detergent wash solution, per detergent manufacturer's instructions, using lukewarm tap water. Submerge and soak instruments in wash solution for a minimum of 1 minute.
 - c. While still submerged, remove visible soil by scrubbing with a soft, nylon bristle brush for a minimum of 4 minutes or until no visible soil is observed.
 - d. Use appropriately sized lumen brush to scrub the interior of lumen and soft bristled, nylon brush to brush the exterior of the lumen. Brush for a minimum of 1 minute, or until no soil is visible.
2. Rinse with running, lukewarm, DI, RO, or filtered water for a minimum of 1 minute for each instrument.
3. Load instruments into washer/disinfector in accordance with the manufacturer's instructions.
 - Arrange instruments with curved surfaces and lumens facing downward to prevent pooling of water on the instrument.
4. Operate the washer/disinfector cycle according to the manufacturer's instructions.
 - Recommended minimal Automated washer/disinfector cycle parameters:
 - a. Pre-wash rinse with cold tap water for 2 minutes.
 - b. Heated wash at 140°F (60°C) for 2 minutes using an enzymatic, pH neutral cleaner such as Prolystica Ultra Concentrate Enzymatic Cleaner.
 - c. Heated tap water rinse at 140°F (60°C) for 20 seconds
 - d. Heated deionized water rinse at 180°F (82°C) for 2 minutes
 - e. Forced air drying at 240°F (116°C) for 9 minutes
 - f. If any residual moisture is observed, wipe dry with absorbent, lint free cloth, or pressurized air for lumens.

9. POST CLEANING INSPECTION:

After cleaning, visually inspect each test article with the naked eye under normal lighting conditions to determine if all visible soil (e.g. blood protein substances and other debris) has been removed. If any soil is still visible, repeat cleaning steps.

9.1 SURGICAL ARM INSPECTION:

Articulating Arms

- Inspect entire assembly for damage.
- Hold arm assembly at column and turn central tightening knob clockwise.
- Check to make sure that arm is rigid at all three joints.
- Insert arm column into table clamp, turn column tightening lever clockwise and ensure that it holds securely.

Note: Loosen articulating arms prior to sterilization.

Flexible Arms

- Inspect entire assembly for damage.
- Turn the flex arm tightening knob clockwise and ensure the arm is sufficiently rigid for intended use.
- When loosened, check cable between links for fraying. Normal use will eventually cause wear to the steel tensioning cable. If cable shows any frayed or broken wires, the flex arm needs to be replaced immediately.

10. STERILIZATION INSTRUCTIONS:

10.1 LUBRICATE:

For instruments with moving parts, lubricate joints with a steam-permeable, water soluble instrument lubricant prior to sterilization.

10.2 STERILIZE:

- Instruments should be sterilized in the open or unlocked position. Central knob of articulating arms must be opened for sterilization.
- Instruments should be sterilized by standard cycles using steam with established procedures.
- We recommend the following sterilization temperature and time.

Gravity

The gravity displacement autoclave process is to sterilize the instruments at 250°F (121°C) for 30 minutes with a 150 minute dry time.

Prevacuum United States Standards

The prevacuum autoclave process is to sterilize the instruments at 270°F (132°C) for 4 minutes with a 30 minute dry time.

Prevacuum EU Standards

The prevacuum autoclave process used to sterilize instruments according to EU standards is at 273°F (134°C) for 18 minutes with a 20 minute dry time.

Autoclave temperatures should not surpass 280°F (137°C), as the handle, insulation or other non-metallic parts may be affected. (Note: The steam autoclave manufacturer may be contacted to confirm appropriate temperatures and sterilization times.)

11. STORE:

Instruments should be stored in a clean dry area with tip protectors.

Please examine instrument prior to use for functionality and damage. When necessary, dispose of products in accordance with national regulations and approved hospital practices for surgical instrumentation disposal.

12. SERIOUS INCIDENTS REPORTING:

Any serious incidents that has occurred in relation to the device should be reported to TSI and the Competent Authority of the Member State in which the user and/or patient is established. Adverse events and incidents can be reported to TSI per contact information on the device labeling.

13. PRODUCT LIFETIME:

End of life is normally determined by wear and damage due to use. Refer to the assembly instructions above to ensure that the products function as outlined.

14. WARNING/PRECAUTIONS:

1. CAUTION: US Federal law restricts this device to sale by or on the order of a physician.
2. Product is intended to be used by trained surgeons.
3. TSI products are to be used only with the TSI retractor systems and may not be used with other manufacturer's products.
4. Items modified or repaired outside of TSI may result in damage or loss of functionality to a device. TSI is then, no longer responsible for failure of the device.
5. Use of this instrument for any purpose, or in any manner other than those described here may cause instrument damage or failure which could result in serious patient injury or death. If needed, all TSI metal products or fragments thereof may be located by means of an X-Ray.
6. To prevent corrosion, instruments made of different alloys should be physically separated during cleaning and sterilization.
7. Products must be in unlocked or loosened position prior to sterilization.
8. For instruments with moving parts, lubricate joints with a steam-permeable, water soluble instrument lubricant prior to sterilization.
9. CJD (Creutzfeldt-Jakob Disease): Discard or destroy any product that comes in contact with or is exposed to patients with CJD, or anyone suspected of CJD. TSI does not provide any validated instructions to eliminate risk of cross-contamination.
10. To maintain intended clamping capacity of the table clamp, do not tighten the rail clamping knob when the articulating arm column is not fully installed.
11. The light source must remain off until the light cable is inserted into the retractor blade(s).
12. Place the light source away from items that are flammable.
13. TSI light cables should only be used with the TSI light source (ML-0051).
14. Reusable Light Cables Only: To verify that the proper amount of light output is achieved, hold single fiber optic end of light cable up to room light and look in bifurcated end to check for the percentage black dots seen (the black dots represent broken fibers in the bundle). If greater than fifty percent (50%) of the fibers are broken, the light cable should be replaced.
15. Once the light cable is connected to the light source, do not place the light cable on drapes, sponges, or any flammable object.
16. Once the light cable is connected to the light source, do not allow the light cable to hang over the side of the sterile field.
17. To prevent damage to the screws and forks on the connection points of the Retractor Ring, assemble and disassemble ring on a flat surface.
18. Adjusting any Flex Arm's position without loosening tension knob will cause cable wear.
19. Acute bending of any Flex Arms will cause cable damage/ wear.
20. DO NOT FORCE ANY KNOBS PAST STOP.
21. Single use devices may not be reused. Reuse may compromise structural integrity

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