ATP STERILE SCREWS INSTRUCTIONS FOR USE
VIS STÉRILES ATP MODE D'EMPLOI
VITI STERILI ATP ISTRUZIONI PER L'USO
ATP STERILE SCHRAUBEN ANWEISUNGEN ZUR VERWEN-
DUNG
ATP STERIELE SCHROEVEN INSTRUCTIES VOOR GEBRUIK
PARAFUSOS ESTÉREIS ATP INSTRUÇÕES DE USO

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UL-0370S-05, UL-0372S-05, UL-0374S-05, UL-376S-05, UL-0378S-05

Company Information:

For more information on this product
or any products mentioned in this IFU
please contact:

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TeDan Surgical Innovations, Inc.

Product Description:

The TeDan Surgical Innovations, Inc. (TSI) Phantom UL ATP Sterile Screws are single use, STERILE devices available in ten screw lengths (60, 80, 100, 120, 140, 160, 180, 200, 220, 240mm). They are supplied sterile, one (1) screw per pouch and 5 pouches per box, and are intended for single use only.

REF#'s:

Product No.	Screw Length (mm)	Qty of Screws/Pouch	Qty of Pouches/Box
UL-03605-05	60	1	5
UL-03625-05	80	1	5
UL-03645-05	100	1	5
UL-03665-05	120	1	5
UL-03685-05	140	1	5
UL-03705-05	160	1	5
UL-03725-05	180	1	5
UL-03745-05	200	1	5
UL-03765-05	220	1	5
UL-03785-05	240	1	5

The TSI Phantom UL ATP Sterile Screws are intended for the fixation of the blade to the vertebral body to prohibit drifting in the Anterior to the Psoas approach to the lumbar spine. Product is intended to be used by trained surgeons and is designed for use with TSI Phantom Screwdrivers and Phantom UL Blade.

Surgical Setting:

Always ensure use of TSI instruments in a strictly sterile clinical environment.

Indication for Use:

Anterior to the Psoas surgery.

Intended Patients Groups:

Patients for which lumbar spine surgery from an anterior to the psoas approach is required.

Intended User Groups:

Licensed practicing physicians who have been informed on the use of Phantom ATP Sterile Screws used with TSI retractor systems. It is also expected that clinical assistants are knowledgeable in setting up the retractor system.

Contraindications:

No known contraindications.

Intended Purpose:

The Phantom ATP Sterile Screws are used in combination with Phantom UL Blades to fixate blades to the vertebral bodies. Products have been designed to be used with Phantom UL Blades and TSI Screwdrivers.

Intended Clinical Benefits:

Provides a temporary anchoring mechanism for the blades during anterior to the psoas lumbar spine surgery.

Instructions for use:

1. Inspect package prior to use and do not use if the package is damaged or sterility has been compromised. If packaging has been damaged or opened before use, DO NOT USE. Report any damaged packaging to TSI.
2. After selecting the desired length of retractor blade, choose the corresponding screw. The length should be selected based upon blade identification.

3. Using a TSI Screwdriver (such as ML-0055), first insert and then screw until seated. The Phantom ATP Sterile Screw should go through a Phantom UL Blade Cannulation.
4. Remove the screwdriver.
5. Repeat steps 1 through 4 if a second Phantom ATP Sterile Screw is desired. Only one screw per cannulated blade is recommend.
6. The above procedure ensures correct placement of the screws. Do not bend the screws as this may cause failure of the screw.
7. Remove the Phantom ATP Sterile Screws after use with a TSI Screwdriver (such as ML-0055).
8. TSI Phantom ATP Sterile Screws are single use only products. After the procedure, dispose of the screws in an approved sharps container.
9. If needed, the screw or any fragment thereof can be located by means of an X Ray.

Reporting Serious Incidents:

Any serious incidents that has occurred in relation to the device should be reported to TSI and the Competent Authority of the Member State in which the user and/or patient is established. Adverse events and incidents can be reported to TSI per contact information on the device packaging.

Warning/Precautions:

1. CAUTION: US Federal law restricts this device to sale by or on the order of a physician.
2. This product is supplied STERILE and Single Use Only. Do not reuse. Reuse may compromise the structural integrity of the device and/or lead to device failure. Reuse also presents biological hazards associated with disease transmission and immune/allergy issues, some of which could cause severe illness or be fatal.
3. Do not bend the screws as this may cause failure of the screw.
4. The decision to use Phantom ATP Sterile Screws is based upon suitability of the patient's bone condition.
5. Phantom ATP Sterile Screws are only compatible with TSI screwdrivers.
6. Do not use this device for any purpose other than the stated intended use.
7. Potential Adverse Events / Residual Risks: Infection / Adverse Reaction, Delay in Surgery, Injury.
8. Use of this instrument for any purpose, or in any manner other than those described here may cause instrument damage or failure which could result in serious patient injury or death. If needed, all TSI metal products or fragments thereof can be located by means of an X-Ray.
9. TSI Phantom ATP Sterile Screws are for use only with the TeDan Phantom UL Retractor system and may not be used with other manufacturer's products.

Sterilization and Shelf Life:

Phantom ATP Sterile Screws are STERILE single use only and are sterilized by ethylene oxide. They may not be re-sterilized or re-used. Sterile packages have a five (5) year shelf life and should be stored at room temperature.

See www.tedansurgical.com/symbolsglossary for the Symbol Glossary.

Company Information:

For more information on this product or any products mentioned in this IFU please contact:


TeDan Surgical Innovations, Inc.
12320 Cardinal Meadow Drive
Suite 150
Sugar Land, TX 77478 USA

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TeDan Surgical Innovations, B.V.